


# MILTEX HI-LIGHT

## STERILE DISPOSABLE BIOPSY PUNCH WITH PLUNGER SYSTEM




The Miltex® Sterile Disposable Biopsy Punch, featuring the patent pending internal plunger system, addresses the challenge of removing a skin specimen when it becomes lodged inside the metal lumen of the punch. The internal plunger system operates similarly to a retractable pen, providing the user with effortless ejection of a lodged skin specimen. This unique feature can help eliminate the problem of "crush" artifact.

Punches are available in five sizes, marked clearly on the handle, blister pack and box. Punches are individually packed in color-coded, transparent, sturdy blister packs, 25 per dispenser box.

REF	Size
33-31AA-P/25	1.0mm w/plunger
33-31A-P/25	1.5mm w/plunger
33-31-P/25	2.0mm w/plunger
33-32-P/25	3.0mm w/plunger
33-34-P/25	4.0mm w/plunger

2 NEW!  
SIZES



*Miltex*  
Redefining Excellence

589 Davies Drive  
York, PA 17402  
phone 717 840-9335  
toll free 866 854-8300  
fax 866 854-8400  
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

0107/5/B/3

SURG-PLUNG-HI




## Sterile Disposable Biopsy Punches

(Without Plunger System)

Now you can have the convenience of disposable punches along with the quality you expect from Miltex, Inc. Miltex offers an extensive product line of sterile disposable biopsy punches that are manufactured from the finest grade of surgical steel to ensure superior cutting, every time.

- Seamless, stainless steel, razor-sharp cutting edges eliminate torn or pulled biopsies
- Ribbed handle provides comfort and control
- Individual color-coded blister packs
- Convenience of disposables, single-use
- Sizes clearly marked on each punch
- Each dispenser box contains 50 punches
- Ten sizes for a wide range of procedures
- Assortment box available and contains 10 each of the following sizes: 2mm, 3mm, 4mm, 5mm, and 6mm



REF	Size
33-31AA	1.0mm
33-31A	1.5mm
33-31	2.0mm
33-31B	2.5mm
33-32	3.0mm
33-33	3.5mm
33-34	4.0mm
33-35	5.0mm
33-36	6.0mm
33-37	8.0mm
33-38	Assorted sizes*

Each box contains 50 sterile color-coded blister packs.

\*10 each of 2mm, 3mm, 4mm, 5mm, and 6mm.