


MILTEX HI-LIGHT

BAGGISH CERVICAL BIOPSY PUNCH FORCEPS



REF	Description
30-1450	BAGGISH Cervical Biopsy Punch Forceps w/lock, 8-1/4" (21 cm), 6 x 4 x 1.5 mm bite

Miltex® Cervical Biopsy Punch Forceps are manufactured to traditional standards of German craftsmanship as well as using only the highest quality stainless steel. Handles are ergonomically designed to provide the clinician with a comfortable grip and positive control, requiring only the slightest pressure to excise cervical tissue. In addition to the Baggish pattern, Miltex offers a variety of instruments, which include different handle configurations, jaw sizes and shapes for all cervical and uterine biopsy needs.

REF	Description
30-1442	TISCHLER Cervical Biopsy Punch Forceps, 9-3/4" (24.8 cm) shaft, 6 X 3 X 1.5 mm bite
30-1443	BURKE (Modified Baby TISCHLER) Cervical Biopsy Punch Forceps, 7-3/4" (19.7 cm) shaft, 4.5 X 3 X 1.5 mm bite
30-1445	TOWNSEND Mini-Bite Cervical Biopsy Punch Forceps, 7-1/2" (19.1 cm) shaft, 4.2 X 2.3 X 1 mm bite, w/lock
3-700	MILTEX Spray Lube, 8 oz. (0.24 liter) bottles (case of 12). Pre-mixed, non-silicone, anti-corrosive formula in a ready-to-use 8 oz. spray pump bottle. Steam permeable and water soluble. Helps prevent spotting, staining and rusting. Lubricates instruments.

It is important to lubricate biopsy punch forceps prior to sterilization.

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

Miltex
Redefining Excellence

www.miltex.com
toll free 866 854-8300
phone 717 840-9335