

MILTEX HI-LIGHT

CERVICAL BIOPSY PUNCH FORCEPS



30-1442



30-1443



30-1445



30-1450

Miltex offers an extensive selection of finely crafted Cervical and Uterine Biopsy Punch Forceps which consistently provide a sharp cut and artifact free biopsy. These instruments are guaranteed to be as durable as they are exact in the results they provide for patient evaluation. Miltex German crafted Gyn Biopsy Punch Forceps exceed standards for providing exceptional trouble-free performance and reliable clinical results.

It is critical that the Clinician have a range of choices to optimize biopsy results and Miltex offers numerous jaw configurations to support this need. Handles are ergonomically designed to provide the clinician with a comfortable grip and positive control, requiring only the slightest pressure to excise cervical or uterine tissue.

REF	Description
30-1442	TISCHLER Cervical Biopsy Punch Forceps, 9-3/4" (24.8 cm) shaft, oblong 7 X 3 X 1.5 mm bite
30-1443	BURKE (Modified Baby TISCHLER) Cervical Biopsy Punch Forceps, 7-3/4" (19.7 cm) shaft, oval 4.5 X 3 X 1.5 mm bite
30-1445	TOWNSEND Mini-Bite Cervical Biopsy Punch Forceps, 7-1/2" (19.1 cm) shaft, 4.2 X 2.3 X 1 mm bite, w/lock
30-1450	BAGGISH Cervical Biopsy Punch Forceps, 8-1/4" (21 cm), oval 6 x 4 x 1.5 mm bite, w/lock

It is important to lubricate biopsy punch forceps prior to sterilization.

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

Miltex
Redefining Excellence


589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

1006/7.5/A/1

SURG-GYNPNCN-HI




MILTEX® UTERINE BIOPSY PUNCH FORCEPS



REF	Description
30-1480	MILTEX Uterine Biopsy Forceps, 9" (22.9 cm) shaft, round cup jaw, 6 mm bite
30-1482	KEVORKIAN-YOUNGE Uterine Biopsy Forceps, 9-1/2" (24.1 cm) shaft, rectangular 3.5 x 8 mm bite, teeth on lower jaw
30-1485	EPPENDORFER Uterine Biopsy Forceps, 8-1/2" (21.6 cm) shaft, 4 x 4.5 x 8 mm bite
30-1486	WITTNER Uterine Biopsy Forceps, 8-1/2" (21.6 cm) shaft, 3 x 4.5 x 8 mm bite with teeth on lower jaw, angled tip
30-1487	WITTNER Uterine Biopsy Forceps, 8-1/2" (21.6 cm) shaft, 3 x 4.5 x 8 mm with teeth on lower jaw, straight tip

MILTEX® ROTATING SHAFT UTERINE BIOPSY FORCEPS

9" (22.9 cm) shafts are interchangeable on Pistol Grip Spring Handles



The Miltex two piece fully-rotational forceps permit improved site visualization for the Clinician and unimpaird access to the target tissue. Interchangeable shafts with different jaw sizes and bite patterns are available to insure that an appropriate amount of tissue is collected while minimizing patient discomfort.

REF	Description
30-1490	KEVORKIAN-YOUNGE Biopsy Forceps, rectangular 3.5 x 8 mm bite, teeth on lower jaw, complete with pistol grip handle
30-1491	KEVORKIAN-YOUNGE Biopsy Forceps, shaft only
30-1492	EPPENDORFER-KRAUSE Biopsy Forceps, oval 4.5 x 6.5 mm bite, complete with pistol grip handle
30-1493	EPPENDORFER-KRAUSE Biopsy Forceps, shaft only
30-1494	TISCHLER Biopsy Forceps, oblong 7 x 3 x 1.5 mm bite, complete with pistol grip handle
30-1495	TISCHLER Biopsy Forceps, shaft only
28-400	Pistol Grip Spring Handle for MILTEX Uterine Biopsy Forceps 30-1490 to 30-1495