


MILTEX HI-LIGHT

COLPOSCOPY SETS



Standard Colposcopy Set

30-15	GRAVES Vaginal Speculum, medium, 1-3/8" (3.5 cm) X 4" (10.2 cm)
30-55	PEDERSON Vaginal Speculum, narrow blades, medium, 7/8" (2.2 cm) X 4" (10.2 cm)
30-951	EMMETT Uterine Tenaculum Hook, 9" (22.9 cm) sharp points, style 1, half-curved, small
30-1350	KOGAN Endocervical Speculum, 9-1/2" (24.1 cm), w/ gauge & set screw, fenestrated jaws 7 X 25 mm
30-1382	KEVORKIAN-YOUNGE Endocervical Biopsy Curette, 12" (30.5 cm), loop 3 X 12 mm, with basket
30-1442	TISCHLER Cervical Biopsy Punch Forceps, 9-3/4" (24.8 cm) shaft, 7 x 3 x 1.5 mm bite

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

See reverse side for additional information


0307/2/C/3

SURG-COLP-SET



+H834SURGCOLPSET0Q

COLPOSCOPY SETS



Delicate Colposcopy Set

30-15	GRAVES Vaginal Speculum, medium, 1-3/8" (3.5 cm) X 4" (10.2 cm)
30-55	PEDERSON Vaginal Speculum, narrow blades, medium, 7/8" (2.2 cm) X 4" (10.2 cm)
30-951	EMMETT Uterine Tenaculum Hook, 9" (22.9 cm) sharp points, style 1, half-curved, small
30-1351	MILTEX-KOGAN Endocervical Speculum, 11" (27.9 cm), ratchet & gauge, extra delicate fenestrated jaws 3 x 25 mm
30-1383	KEVORKIAN-YOUNGE Endocervical Biopsy Curette, 12" (30.5 cm), loop 3 X 12 mm, without basket
30-1443	BURKE (BABY-TISCHLER) Cervical Biopsy Punch Forceps, 7-3/4" (19.7 cm) shaft, 4.5 x 3 x 1.5 mm bite

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex®, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.