

# MILTEX® HI-LIGHT

## SURGICAL INSTRUMENT CLEANER & LUBRICANT



3-710

### INSTRUMENT CLEANER & LUBRICANT

A dual-purpose instrument cleaner and lubricant. Chemically engineered for use as a presoak or for manual cleaning, washer/decontaminators, ultrasonic cleaners and cart washers.

- A special penetrating lubricant is blended into the cleaning formula, protecting the instruments and extending their life.
- Detergents and emulsifiers quickly break up and remove organic soil from box locks and other hard to reach places restoring the look and feel of a new instrument.

REF	Description
3-710	1 Gallon (128 fl. oz.) 3.79L (see label for concentration rates)

*Miltex*  
Redefining Excellence

589 Davies Drive  
York, PA 17402  
phone 717 840-9335  
toll free 866 854-8300  
fax 866 854-8400  
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email [customerservice@miltex.com](mailto:customerservice@miltex.com), or visit us at [www.miltex.com](http://www.miltex.com).

0507/5/A/1

SURG-CLNLUB-HI



\*+H834SURGCLNLUBHI0Q\*

# INSTRUMENT LUBRICANT



AVAILABLE  
SUMMER '07

3-705

## INSTRUMENT LUBRICANT

A water emulsible soaking agent which works to extend the life of surgical instruments. Concentrated non-silicone lubricant and rust inhibitor. Used daily and prior to sterilization to keep moving parts operating smoothly.

- Formulated to extend the working life of surgical instruments.
- Helps Prevent Corrosion.
- Formulated with gentle, non-aggressive ingredients.
- Frees Up Box Locks – Lubricates moving parts of surgical instruments.

REF	Description
3-705	1 Gallon (128 fl. oz.) 3.79L (see label for concentration rates)