


MILTEX® HI-LIGHT

NEW AND IMPROVED McGIVNEY HEMORRHOIDAL LIGATOR

(FOR USE WITH BOTH LATEX-FREE & STANDARD LATEX O-RINGS)



28-154B (includes loading cone REF# 28-156)

28-155

28-152

PRODUCT DESCRIPTION

Improved lightweight model with comfortable offset handle for controlled grip and better viewing access. The shaft with ligating drum can be rotated 360 degrees for correct positioning and is easy to disassemble for thorough cleaning. Stainless Steel loading cone facilitates placement of the O-ring bands onto the ligating drum. Compatible with Latex and Latex-Free ligating O-rings. (See reverse side for additional product details.)

REF

28-154A	McGIVNEY Hemorrhoidal Ligator Set for LATEX-FREE applications - includes McGivney Hemorrhoidal Ligator 8" shaft (PN 28-154B), McGIVNEY Grasping Forceps 7-1/2" (PN 28-152), Loading Cone (PN 28-156) and 100 LATEX-FREE O-Rings (28-155).
28-154B	McGIVNEY Hemorrhoidal Ligator improved model for LATEX-FREE applications - 8" (20.3 cm) working length. Supplied with Loading Cone (28-156) and 100 LATEX-FREE O-Rings (28-155).
28-152	McGIVNEY Hemorrhoidal Grasping Forceps, 7-1/2" (19.1 cm), angled shanks. Specially designed for use with McGivney Ligators.
28-155	LATEX-FREE O-Rings compatible with PN 28-154B Ligator only. 100 rings per bag. Provided Non-Sterile.
28-151	LATEX O-Rings compatible with PN's 28-150, 28-150XL & 28-154B Ligator patterns. 100 rings per bag. Provided Non-Sterile. (See reverse side for additional product details.)
28-156	Loading Cone compatible with PN 28-154B Ligator only.

See reverse side for additional product details

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.



589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

