



# MILTEX HI-LIGHT


## MICRO ALLIGATOR EAR FORCEPS




**19-2081** (Satin Finish) **NEW**  
**19-2081-B** (Ebony Finish)  
Micro Alligator Ear Forceps  
3-1/4" (8.3cm) Shaft  
0.6mm wide Serrated Jaws




**19-2090** (Satin Finish) **NEW**  
**19-2090-B** (Ebony Finish)  
Micro Alligator Ear Forceps  
3-1/4" (8.3cm) Shaft  
0.5mm wide Oval Cup Jaws



**19-2100** (Satin Finish)  
**19-2100-B** (Ebony Finish)  
Micro Alligator Ear Forceps  
3-1/4" (8.3cm) Shaft  
0.8mm wide Serrated Jaws



**19-2110** (Satin Finish)  
**19-2110-B** (Ebony Finish)  
Micro Alligator Ear Forceps  
3-1/4" (8.3cm) Shaft  
0.8mm wide Oval Cup Jaws



### FEATURES

- The most delicate instruments for otological surgery
- Excellent for foreign body removal
- Ideal for precision grasping and extraction of cerumen from the ear canal
- Serrated and Oval cup jaw patterns available (0.5mm to 0.8mm wide)
- Available in either ebonized or stainless steel finish - Autoclavable
- Manufactured to traditionally high Miltex standards and German craftsmanship

*Miltex*  
Redefining Excellence

589 Davies Drive  
York, PA 17402  
phone 717 840-9335  
toll free 866 854-8300  
fax 866 854-8400  
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email [customerservice@miltex.com](mailto:customerservice@miltex.com), or visit us at [www.miltex.com](http://www.miltex.com).

0404/2/B/1

SURG-EB-EAR



\*+H834SURGEBEARZ\*