

NEW PRODUCT HI-LIGHT

MINI-OSTEOTOMES

Ideal for maxillofacial, oral, plastic and reconstructive surgery procedures



REF	Description	Length	Width
PM-044402	Mini Osteotome, Straight	5" (127mm)	2mm
PM-044404	Mini Osteotome, Straight	5" (127mm)	4mm
PM-044406	Mini Osteotome, Straight	5" (127mm)	6mm
PM-044408	Mini Osteotome, Straight	5" (127mm)	8mm
PM-044410	Mini Osteotome, Straight	5" (127mm)	10mm
PM-044412	Mini Osteotome, Straight	5" (127mm)	12mm
PM-044422	Mini Osteotome, Curved	5" (127mm)	2mm
PM-044424	Mini Osteotome, Curved	5" (127mm)	4mm
PM-044426	Mini Osteotome, Curved	5" (127mm)	6mm
PM-044428	Mini Osteotome, Curved	5" (127mm)	8mm
PM-044430	Mini Osteotome, Curved	5" (127mm)	10mm
PM-044432	Mini Osteotome, Curved	5" (127mm)	12mm

Distal working end serrated for bone compacting and mallet interface control

EAR LOBE FORCEPS

Ideal for the repair of torn or stretched ear lobe piercings



REF	Description	Length
PM18-1220	Ear Lobe Forceps, with set screw, 8 x 10 mm inside ring diameter	3-1/2" (8.9 cm)

See Reverse Side For Additional New Padgett Instruments

Illustrations and content signify general description only and may be subject to minor changes. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 1-866-854-8300, e-mail customerservice@miltex.com, or visit us at www.miltex.com.



+H834SURGNEWPADGH10+

Miltex
Redefining Excellence

589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

SKIN GRAFT KNIVES



PM-14701

PM-4597

PM-4597LH

PM-4598

PM-4596

PM-8632

The Silver Miniature Skin Graft Knife is ideal for office use involving small "postage stamp" size grafts. It can also be used for wound debridement.

The Watson Knife is a modified Braithwaite design with a fine smooth rod which glides on the skin. It is a versatile and accurate hand held knife with an adjustable cut depth of 0.1mm to 1.0mm.

The Cobbett Knife is also a modified Braithwaite design with a grooved rod on which the knife oscillates.

USF Debriding Forceps have angled jaws for optimized visualization and control during procedures. They are specially designed for blunt debridement of a burn or wound. The blunt tip and side blades help prevent damage to viable skin.

REF	Description	Length	Blade
PM-14701	Silver Miniature Skin Graft Knife	7-1/2" (19.1cm)	uses standard double-edge razor blade
PM-4597	Watson Skin Graft Knife*	12" (30.5cm)	uses PM-4596 blades
PM-4597LH	Watson Skin Graft Knife, Left-Hand*	11-3/4" (29.8cm)	uses PM-4596 blades
PM-4598	Cobbett Graft Knife*	12-1/2" (31.8cm)	uses PM-4596 blades
PM-4596	Sterile Blades, 10 pack	6-1/4" (15.9cm)	for skin graft knives PM-4597, PM-4597LH, PM-4598
PM-8632	USF Debriding Forceps, Angled Jaws	7" (17.8cm)	

* Modified Braithwaite design