

MILTEX® NEW PRODUCT HI-LIGHT

WING-GRIP™ NAIL NIPPERS WITH AN INTERNAL PIN MECHANISM!

Precision Jaw
Exact blade alignment
for effortless cutting.

Internal Pin Mechanism
enhances the performance
of the instrument.

"Wing" Design Handle
comfortable "grip" and
optimal performance.



The new Miltex® Wing-Grip™ nippers feature a jaw design with exact blade alignment which facilitates effortless cutting and extends the performance life of the instrument. The unique ergonomic, wing design handle combines shape, length, weight, balance and feel for optimal comfort and performance. A fully tempered box lock maintains precise instrument alignment. The new Miltex® nippers feature an internal pin mechanism. Two sizes are available.

REF	Description
40204	Wing-Grip™ Nail Nipper 4-1/2" (11.4 cm), curved jaws
40208	Wing-Grip™ Nail Nipper 5" (12.7 cm), curved jaws

Illustrations and content provide general description only and may be subject to change. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 1-866-854-8300, e-mail customerservice@miltex.com, or visit us at www.miltex.com.

Miltex and Miltex Redefining Excellence are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation.
Wing-Grip is a trademark of Miltex, Inc. a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra Life Sciences Corporation.
Copyright © 2010 Integra LifeSciences Corporation

