


# MILTEX HI-LIGHT



## ROOT PLANE PETITE CURETTES



The Miltex GripLite™ Root Plane Petite Curettes are area specific instruments designed after standard Gracey instruments. The terminal shank is 3mm longer which allows for access into deeper pockets and root surfaces. The blade is half the length of a standard Gracey instrument and root plane instruments, allowing for better scaling in narrow pockets and around tooth furcations.

| REF     | Description                              |
|---------|--|
| 70-5119 | 1/2 Root Plane Petite GripLite Curette   |
| 70-5120 | 3/4 Root Plane Petite GripLite Curette   |
| 70-5121 | 5/6 Root Plane Petite GripLite Curette   |
| 70-5122 | 7/8 Root Plane Petite GripLite Curette   |
| 70-5123 | 9/10 Root Plane Petite GripLite Curette  |
| 70-5124 | 11/12 Root Plane Petite GripLite Curette |
| 70-5125 | 13/14 Root Plane Petite GripLite Curette |
| 70-5126 | 15/16 Root Plane Petite GripLite Curette |
| 70-5127 | 17/18 Root Plane Petite GripLite Curette |
| 70-5500 | Color Code Rings, Kit                    |



70-5122



70-5123



70-5124



70-5119



70-5125



70-5120



70-5126



70-5121



70-5127



www.miltex.com  
toll free 866 854-8300  
phone 717 840-9335

All instruments available for individual sale. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email [customerservice@miltex.com](mailto:customerservice@miltex.com), or visit us at [www.miltex.com](http://www.miltex.com).

