

MILTEX® HI-LIGHT

GRIP DEEP™ EXTRACTING FORCEPS VETERINARY DENTISTRY



Miltex® Grip Deep™ Extracting Forceps are designed to allow subgingival access to the root during extractions. Tapered beaks can provide deeper access into the tooth socket and maximize the tooth-to-forceps contact. Serrated beaks aid in gripping the tooth surface which may assist in atraumatic removal of the tooth and can potentially reduce the risk of fractures, or causing trauma to the surrounding hard and soft tissues.

The take-apart feature of the Miltex® Grip Deep™ Extracting Forceps have a right and left side which can be separated to allow access between the box-lock, hinge area for instrument cleaning and processing.

Available in 11 patterns.

REF	Description	Length
DEFDG1STD	Extracting Forceps, Grip Deep 1 Standard	6-1/4"
DEFDG99C	Extracting Forceps, Grip Deep 99C	6-1/4"
DEFDG150	Extracting Forceps, Grip Deep 150	6-1/4"
DEFDG150A	Extracting Forceps, Grip Deep 150A	6-1/4"
DEFDG151	Extracting Forceps, Grip Deep 151	6-1/4"
DEFDGM3	Extracting Forceps, Grip Deep MD3	6-1/4"
DEFDGM4	Extracting Forceps, Grip Deep MD4	6-1/4"
DEFDG210S	Extracting Forceps, Grip Deep 210S	6-1/4"
DEFDG222	Extracting Forceps, Grip Deep 222	6-1/4"
DEFDG65	Extracting Forceps, Grip Deep 65	6-1/4"
DEFDG845	Extracting Forceps, Grip Deep 845	6-1/4"



Take-apart feature!

See reverse side for beak details

Illustrations and content signify general description only and may be subject to minor changes. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 866-854-8300, e-mail: customerservice@miltex.com, or visit us at www.miltex.com.


Miltex is a registered trademark of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. Miltex Redefining Excellence and Grip Deep are trademarks of Miltex, Inc. The Integra Wave Logo is a trademark of Integra Life Sciences Corporation. Copyright © 2008 Integra LifeSciences Corporation



+H834VETGRPDPHI0U

Miltex
Redefining Excellence
INTEGRA™

toll free 866 854-8300
fax 866 854-8400
www.miltex.com



Grip Deep™ Extracting Forceps 150 shown.

Clinical photography provided by: *Barden A. Greenfield, DVM, Southwind Animal Hospital, PLC*
www.southwindvets.com



DEFDG1STD
Grip Deep 1 Standard



DEFDGM4
Grip Deep MD4



DEFDG99C
Grip Deep 99C



DEFDG210S
Grip Deep 210S




DEFDG150*
Grip Deep 150
** Most popular pattern*



DEFDG222
Grip Deep 222



DEFDG150A
Grip Deep 150A



DEFDG65
Grip Deep 65
No serrations



DEFDG151
Grip Deep 151



DEFDGM3
Grip Deep MD3



DEFDG845
Grip Deep 845
No serrations