

MILTEX HI-LIGHT

VETERINARY INSTRUMENT PACKS

FOR VETERINARY DENTAL AND SURGICAL PROCEDURES

WELDIN VETERINARY PERIODONTAL INSTRUMENT KIT - REF# 6809

Kit includes 1 each of:

REF	Description
SVET-5450	GripLite™ S6 Mirror Handle, single-ended, cone-socket
67-696/3*	#3 Cone Socket Front Surface Mirror, 1 mirror
67-696/4*	#4 Cone Socket Front Surface Mirror, 1 mirror
SVET-23E/18P	GripLite™ S6 #23 Explorer/18mm Probe, double-ended
SVET-5004	GripLite™ S6 J33/H5 (Jacquette J33/1YS Scaler)
SVET-767	GripLite™ S6 Morse 0-00 Scaler
SVET-5201	GripLite™ S6 Columbia 13/14 Universal Curette
SVET-5202	GripLite™ S6 Barnhart 1/2 Universal Curette
DTST**	Delron Sharpening Test Sticks, Clear, 1 stick
STN-ARK-CON	Conical Arkansas Sharpening Stone, 1 stone
STN-ARK-FLT-4	#4 Flat Arkansas Sharpening Stone, 1 stone
3-084107	Miltex® Sterilization & Storage Cassette (8" L x 4-7/8" W x 1-1/8" H)

*(not sold individually – REF is 12/box)

** (not sold individually – REF is 6/pk)

Check with distributor for individual sale.



Includes Training DVD & Instructional Chart!

NIEMIEC VETERINARY EXTRACTION INSTRUMENT KIT - REF# 6806

Kit includes 1 each of:

REF	Descriptions
4-7	# 3 Scalpel Handle
73-64	BABY MOLT Periosteal Elevator
DEFF	Extracting Forceps, Universal, Mini
5D-324	LA GRANGE Scissors 4-1/4" Curved, One Serrated Blade, Angled Shank
8-15TC	OLSEN-HEGAR Needle Holder, Combined Suture Scissors 5-1/2"
6-124	BROWN-ADSON Tissue Forceps, 4-3/4", 7x7 Teeth
5-16	Operating Scissors (suture scs), 5-1/2" Straight, S/B
DELLUXS2	Miltex Luxating Elevator, Straight, 2mm
DELLUXS3	Miltex Luxating Elevator, Straight, 3mm
DELLUXS5	Miltex Luxating Elevator, Straight, 5mm
DELLUXC3	Miltex Luxating Elevator, Curved, 3mm
DELLUXC5	Miltex Luxating Elevator, Curved, 5mm
3-085009	Miltex #9 Sterilization & Storage Cassette, 8" L x 6" W x 1-3/16" H



Includes Training DVD & Instructional Guide!

CANINE / FELINE SPAY PACK

Each set includes:

Qty.	REF	Description
1	4-7	#3 Scalpel Handle
1	5-124	Mayo Scissors 6 3/4" Straight
1	6008	Snook Ovariectomy Hook
2	7-38	Kelly Forceps 5 1/2" Curved
1	5-182	Metzenbaum Scissors 7" Curved
1	6-44	Tissue Forceps 5 1/2" 1x2 Teeth
1	16-6	Allis Tissue Forceps 6" 4x5 Teeth
4	7-504	Backhaus Towel Clamps 3 1/2"
1	6-124	Brown-Adson Forceps 4 3/4" 7x7 Teeth
1	8-16TC	Olsen-Hegar Needle Holder 6 1/2" (canine)
or		
	8-15TC	Olsen-Hegar Needle Holder 5 1/2" (feline)
2	7-4	Halsted Mosquito Forceps 5" Curved
2	7-172	Roch-Carmalt Forceps 6 1/4" Curved
1	3-6815	Sterilization/Storage Cassette



15" x 8 1/4" x 1 1/4"

	REF#	REF#
	Canine	Feline
Miltex®	6810	6820
MeisterHand™	MH96810	MH96820
Vantage®	V96810	V96820



589 Davies Drive
York, PA 17402
phone 717 840-9335
toll free 866 854-8300
fax 866 854-8400
www.miltex.com

All instruments available for individual sale. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

0507/5/C/3

VET-PACK-HI



+H834VETPACKHI0Y

MILTEX WINGED ELEVATOR BASIC KIT

Kit is ideal for practices with companion-animal patients.

REF# DELWBASIC - Miltex Winged Dental Elevator Basic Kit Contains:


DELWNG2	1 – Miltex Winged Dental Elevator, Size 2mm
DELWNG3	1 – Miltex Winged Dental Elevator, Size 3mm
DELWNG4	1 – Miltex Winged Dental Elevator, Size 4mm
DELWNG5	1 – Miltex Winged Dental Elevator, Size 5mm
3-083005	1 – Miltex #5 Sterilization & Storage Cassette (8" x 3-13/16" x 1-1/8")



GENERAL SURGERY PACK

Set includes:

Qty.	REF	Description
1	5-124	Mayo Scissors 6 3/4" Straight
1	4-7	#3 Scalpel Handle
1	6008	Snook Ovariectomy Hook
3	7-4	Mosquito Hemostat Forceps 5" Curved
2	7-36	Kelly Hemostat 5 1/2" Straight
1	5-182	Metzenbaum Scissors 7" Curved
1	6-124	Brown-Adson Forceps 4 3/4" 7x7 Teeth
4	7-172	Rochester-Carmalt Forceps 6 1/4" Curved
4	7-506	Backhaus Towel Clamps 5 1/4"
1	8-15TC	Olsen-Hegar Needle Holder 5 1/2"
1	3-6815	Sterilization/Storage Cassette




15" x 8 1/4" x 1 1/4"

	REF#
Miltex®	6800
MeisterHand™	MH96800
Vantage®	V96800

ORTHOPEDIC PACK

Set includes:

Qty.	REF	Description	REF#
1	9-40-20	Suture Wire, 20 Gauge	Miltex® 6840
1	9-40-30	Suture Wire, 30 Gauge	
1	9-124	Wire Cutter	
1	11-74	Senn Retractor, Sharp	
1	11-76	Senn Retractor, Blunt	
1	11-608	Weitlaner Retractor	
1	18-1968	Freer Elevator	
1	19-708	Bone Curette	
1	19-792	Mallet	
1	21-369	Bone File	
1	22-500TC	Wire Twister	
2	26-120	Gigli Handles	
1	26-130	Gigli Saw	
1	27-4	Kern Bone Holding Fcps, Lg	
1	27-7	Kern Bone Holding Fcps, Sm	
1	27-158	Pin Cutter	
1	27-1390	Bone Cutting Forceps	
1	27-603	Osteotome	
1	3-6845	Sterilization/Storage Cassette	



10 3/4" x 8 1/4" x 2 3/4"

MILTEX WINGED ELEVATOR DELUXE KIT

Kit is ideal for practices with diverse patients; companion- to large-animal patients.

REF# DELWDELUX - Miltex Winged Dental Elevator Deluxe Kit Contains:

DELWNG2	1 – Miltex Winged Dental Elevator, Size 2mm
DELWNG3	1 – Miltex Winged Dental Elevator, Size 3mm
DELWNG4	1 – Miltex Winged Dental Elevator, Size 4mm
DELWNG5	1 – Miltex Winged Dental Elevator, Size 5mm
DELWNG6	1 – Miltex Winged Dental Elevator, Size 6mm
DELWNG7	1 – Miltex Winged Dental Elevator, Size 7mm
DELWNG8	1 – Miltex Winged Dental Elevator, Size 8mm
3-084007	1 – Miltex #7 Sterilization & Storage Cassette (8" x 4-7/8" x 1-1/8")




EYE PACK

REF#
Miltex® 6830

Set includes:

Qty.	REF	Description
1	4-7	#3 Scalpel Handle
1	18-35	Wire Eye Speculum
1	18-780	Eye Dressing Forceps 4"
1	18-1462	Stevens Tenotomy Scissors 4 1/8" Straight
1	18-1202	Chalazion Forceps
1	18-1828	Castroviejo Needle Holder 5 1/2" Straight
1	3-6835	Sterilization/Storage Cassette




6 1/2" x 9 1/4" x 1 1/4"

MILTEX® STERILIZATION & STORAGE CASSETTES

Universal Cassette without Mat

Item#	Size
4-6835	6-1/2 x 6-1/4 x 1-3/16
4-083000	8 x 3-13/16 x 1-3/16
4-084000	8 x 4-7/8 x 1-3/16
4-085000	8 x 6 x 1-3/16
4-089100	8 x 9 x 1-3/16
4-008122	8 x 11 x 1-3/16
4-008115	8 x 11 x 2
4-009126	9-1/8 x 13 x 1-3/16
4-6815	15 x 8-1/4 x 1-3/16



Silicone Mat Only (Latex Free)

Item#	Size	Description
4-6835MAT*	6-1/2 x 6-1/4	for Cassette 4-6835
4-083000MAT	8 x 3-13/16	for Cassette 4-083000
4-084000MAT	8 x 4-7/8	for Cassette 4-084000
4-085000MAT	8 x 6	for Cassette 4-085000
4-089100MAT*	8 x 9	for Cassette 4-089100
4-008122MAT*	8 x 11	for Cassette 4-008122 & 4-008115
4-009126MAT*	9-1/8 x 13	for Cassette 4-009126
4-6815MAT	15 x 8-1/4	for Cassette 4-6815

*Item coming soon



3-072007

3-083005



REF	Description	Dimensions Inches	Millimeters
3-072007	Slimline 7*	7 x 2 5/8 x 11/16	178 x 67 x 17
3-072014	Slimline 14*	7 x 4 5/8 x 11/16	178 x 115 x 17
3-083005	5 Utility*	8 x 3 13/16 x 1 1/8	204 x 81 x 30
3-084007	7 Utility*	8 x 4 7/8 x 1 1/8	204 x 124 x 28.5
3-085009	9 Utility*	8 x 6 x 13/16	204 x 154 x 30

Other sizes and accessories available - Refer to the Miltex Veterinary, Surgical or Dental catalogs for additional instruments.

*Hand instruments, scissors and rack system not included