

# MILTEX HI-LIGHT

## PROFESSIONAL SPLINTER REMOVAL KIT

IDEAL FOR OFFICE, HOME, TRAVEL  
AND ALL FIRST AID KITS



V96-410

The Liberator is used to approach the embedded fragment by elevating the epidermis. After the fragment is exposed, the Walter Forceps is used for quick and easy splinter removal.

### FEATURES

- Kit comes with a convenient carrying/storage case
- Liberator Tip stores safely in reversible handle when not in use
- Walter Splinter Forceps has a comfortable wide platform with finger grip for firm and stable control
- All stainless steel - autoclavable
- Kit includes one Walter Splinter Forceps 4-1/2" Straight (V96-312) and one Splinter Liberator 5" with reversible handle (V96-400) (can be ordered separately)

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email [customerservice@miltex.com](mailto:customerservice@miltex.com), or visit us at [www.miltex.com](http://www.miltex.com).

*Miltex*  
Redefining Excellence

589 Davies Drive  
York, PA 17402  
phone 717 840-9335  
toll free 866 854-8300  
fax 866 854-8400  
[www.miltex.com](http://www.miltex.com)

0103/2/B/1

SURG-SPLNT-KT



\*+H834SURGSPLNKT5\*