

MILTEX® HI-LIGHT

LED LOUPES LIGHT SOURCE

The Miltex® portable dental and surgical LED light source produces amazing brightness in a completely portable light. The small power pack is lightweight and uses Lithium Ion cells that can be quickly and easily recharged. The portable LED can be attached to most loupes systems.

Features	Benefits
High intensity portable dental and surgical headlight (80,000 Lux @ 14 inches)	Provides pure white light for true color rendition
Innovative power pack design with external battery cells for quick changing State-of-the-art Li Ion batteries (1 included)	Provides 4.5 hours of use on a single charge
Variable light intensity control	Allows for professional to adjust light as needed
Portable Headlamp Weighs only 0.7 ounces	Lightweight reduced pressure
Quick charger - 90% charge in 90 minutes	Reduced charging time
Flip up dental composite filter available	Prevents composite materials from curing



Choose your own
manufacturer
specific mount!
(see reverse
for details)



Each unit comes complete with
1 each of the following:

1. LED Headlight
2. Battery
3. Battery Charger
4. Battery Power Pack (attaches to belt)
5. AC Adapter
6. Attach mount (varies per REF)

See reverse for more information

Illustrations and content provide general description only and may be subject to change. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 1-866-854-8300, e-mail customerservice@miltex.com, or visit us at www.miltex.com.

Miltex and Miltex Redefining Excellence are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra Life Sciences Corporation. All other trademarks are property of their respective owners. Copyright © 2010 Integra LifeSciences Corporation



+H834LPSLGHTHI07

THE MILTEX® LED LOUPES LIGHT AFFIXES EASILY TO YOUR MILTEX® LOUPES!

CHOOSE YOUR OWN MOUNT!

Manufacturer specific mounts come pre-assembled on the lamp head

Use the simple chart below to choose the mount that matches your loupes!

REF	Description	Mount Type
15050	LED Loupes Light Unit	Miltex® loupes mount
15051	LED Loupes Light Unit	Orascoptix® loupes mount
15052	LED Loupes Light Unit	Sheervision® loupes mount
15053	LED Loupes Light Unit	Surgitel® loupes mount
15054	LED Loupes Light Unit	Large Universal mount
15055	LED Loupes Light Unit	Heine® loupes mount

NOTE: Not all mounts match all loupes. Contact your Miltex Representative to ensure compatibility.

THE FOLLOWING REPLACEMENT PARTS ARE ALSO AVAILABLE:



REF 15061
Flip-Up Dental
Composite Filter



REF 15057
Ion Battery



REF 15056
Replaceable Cable

ASSEMBLY IS EASY!



1. Pull off flip up paddle located on the locking knob.



2. Light will mount onto the loupes as the paddle does; simply push the grooved portion in between the locking knob.



3. Clip battery pack onto your belt.



4. Turn on and use!