

MILTEX HI-LIGHT

MX ESONIC

COMPLETE WITH THREE TIPS
UNIVERSAL, PERIODONTAL & VERSATIP™

Tips Deliver Water from Dental Lines to Scaling

Helps flush away debris and calculus from pocket.

Lightweight & Well Balanced

Helps reduce hand and wrist fatigue.

Quick Disconnect Swivel Coupler


Easy on/off connection helps reduce hand and wrist fatigue.



70-2000

Universal and Periodontal Tips

Allows thorough scaling for natural dentition.



VersaTip™

Allows for periodontal work to be completed around implants and composite restorations.



The Miltex® Mx Scaler is a unique lightweight air driven scaler, which can be used to scale calculus and debris on natural dentition, implants and composites. The Mx Scaler is well balanced and features a quick disconnect swivel coupler. It adapts to any standard 4-hole dental hose and is fully autoclavable.

Miltex
Redefining Excellence

www.miltex.com
toll free 866 854-8300
phone 717 840-9335

For additional information, an in-office demonstration, or to place an order, please contact your local Miltext, Inc. authorized Distributor or Miltext, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.