

MILTEX HI-LIGHT

SILITONE COLOR BANDS



- Color bands are expandable from 4mm to 12mm.
- Withstands temperatures up to 450° F.
- Medical grade silicone
- Exercise caution when applying to instruments with sharp points.

Ref	Description
VSTSA	Grey - Pack of 50 bands
VSTSB	Blue - Pack of 50 bands
VSTSD	Gold - Pack of 50 bands
VSTSG	Green - Pack of 50 bands
VSTSK	Black - Pack of 50 bands
VSTSL	Lt Blue - Pack of 50 bands
VSTSN	Brown - Pack of 50 bands
VSTSP	Pink - Pack of 50 bands
VSTR	Red - Pack of 50 bands
VTSU	Purple - Pack of 50 bands
VTSW	White - Pack of 50 bands
VTSY	Yellow - Pack of 50 bands
VST	Kit of 12 colors - 480 pieces
VSTA	Kit of 12 colors - 120 pieces (10 of each color)

Miltex
Redefining Excellence

www.miltex.com
toll free 866 854-8300
phone 717 840-9335

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.