


MILTEX® HI-LIGHT

VETERINARY WINGED DENTAL ELEVATORS



MILTEX® Winged Dental Elevators, developed and approved by Dr. Robert B. Wiggs, D.V.M., Diplomat, American Veterinary Dental College, ease the extraction process and provide the Veterinarian Professional with instrumentation that meets diverse clinical needs.

Elevators are lightweight and ergonomically designed with small handles to improve gripping power, while also allowing increased control and comfort during extraction procedures. Winged Elevators are shaped to improve instrument adaptation to the root structure. When the instruments are used in a side-to-side rocking/twisting motion, blades provide additional leverage, especially when used between teeth in the interdental space. The winged-blades will hook under the buccal bulge, improving the angle of direction in order to lift the tooth and root.

Available in seven sizes and two kit configurations for simplified ordering. For added convenience, kits are pre-packaged in MILTEX® Storage Cassettes.

See reverse side for details and part numbers.

Illustrations and content provide general description only and may be subject to change. For additional information, to schedule an in-office demonstration, or to place an order, please contact your local Miltex authorized Distributor or Miltex Customer Service at 866-854-8300, e-mail: customerservice@miltex.com or visit us at www.miltex.com.

Miltex is a registered trademark and Miltex Redefining Excellence is a trademark of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra Life Sciences Corporation. Copyright © 2008 Integra LifeSciences Corporation

Miltex
Redefining Excellence

INTEGRA


www.miltex.com
toll free 866 854-8300
phone 717 840-9335



MILTEX® Winged Dental Elevators

REF	Description
DELWNG2	MILTEX® Winged Dental Elevator, Size 2 mm
DELWNG3	MILTEX® Winged Dental Elevator, Size 3 mm
DELWNG4	MILTEX® Winged Dental Elevator, Size 4 mm
DELWNG5	MILTEX® Winged Dental Elevator, Size 5 mm
DELWNG6	MILTEX® Winged Dental Elevator, Size 6 mm
DELWNG7	MILTEX® Winged Dental Elevator, Size 7 mm
DELWNG8	MILTEX® Winged Dental Elevator, Size 8 mm

Elevators are available individually



Images shown actual size.

MILTEX® Winged Elevator Basic Kit

Kit is ideal for practices with companion animals.

REF	DELWBASIC	MILTEX® Winged Dental Elevator Basic Kit Contains:
DELWNG2		1 – MILTEX® Winged Dental Elevator, Size 2 mm
DELWNG3		1 – MILTEX® Winged Dental Elevator, Size 3 mm
DELWNG4		1 – MILTEX® Winged Dental Elevator, Size 4 mm
DELWNG5		1 – MILTEX® Winged Dental Elevator, Size 5 mm
3-083005		1 – MILTEX® #5 Storage Cassette 8" L x 3-13/16" W x 1-1/8" H



MILTEX® Winged Elevator Deluxe Kit

Kit is ideal for practices with diverse patients; companion to large animals.

REF	DELWDELUX	MILTEX® Winged Dental Elevator Deluxe Kit Contains:
DELWNG2		1 – MILTEX® Winged Dental Elevator, Size 2 mm
DELWNG3		1 – MILTEX® Winged Dental Elevator, Size 3 mm
DELWNG4		1 – MILTEX® Winged Dental Elevator, Size 4 mm
DELWNG5		1 – MILTEX® Winged Dental Elevator, Size 5 mm
DELWNG6		1 – MILTEX® Winged Dental Elevator, Size 6 mm
DELWNG7		1 – MILTEX® Winged Dental Elevator, Size 7 mm
DELWNG8		1 – MILTEX® Winged Dental Elevator, Size 8 mm
3-084007		1 – MILTEX® #7 Storage Cassette 8" L x 4-7/8" W x 1-1/8" H

