

MILTEX® HI-LIGHT

METRIC SYRINGES



MILTEX® N-TRALIG® SYRINGE FEATURES:

- *Precise 0.2cc of anesthetic delivered per ratchet*
- *Positive ratchet action helps protect against cartridge breakage*
- *Removable hub allows use of the 30° adapter*

Miltex® N-Tralig® Intraligamentary Syringe is handcrafted from high quality German stainless steel, which gives the syringe a greater life expectancy and consistent and precise delivery of 0.2cc of anesthetic.

REF	Description
76-55M	N-Tralig® Intraligamentary Syringe
76-56M	30° Adapter N-Tralig® Barrel Cap

See reverse side for additional patterns



METRIC SYRINGES

ASPIRATING SYRINGES FEATURE:

- Variety of product options – CW and A type syringes
- GripRite® CW Syringes feature:
 - Silicone covered handles and wings
 - Available in blue and red



76-80M

76-30M



76-32M

76-70M

REF	Description
76-80M	Aspirating Syringe, Type "A", 1.8cc
76-30M	GripRite® Aspirating Syringe, Blue Silicone Grips, CW Type, 1.8cc
76-32M	GripRite® Aspirating Syringe, Red Silicone Grips, CW Type, 1.8cc
76-70M	Aspirating Syringe, CW Type, 1.8cc

PETITE ASPIRATING SYRINGES FEATURE:

- Reduced oval shaped thumb ring diameter
- Silicone covered handles and wings
- Flared base on ring handle to ease pressure while dispensing
- Shortened harpoon rod to reduce strain during aspirations



76-42M



76-40M

REF	Description
76-40M	GripRite® Petite Aspirating Syringe, Blue Silicone Grips, CW Type, 1.8cc
76-42M	GripRite® Petite Aspirating Syringe, Red Silicone Grips, CW Type, 1.8cc

SELF ASPIRATING SYRINGE FEATURES:

- Lightweight aluminum handle design
- Blunt harpoon rod which allows aspiration by means of counter pressure



76-90M

REF	Description
76-90M	Lightweight Self-Aspirating Syringe, 1.8cc

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex Customer Service at 1-717-757-7972, email internationalsales@miltex.com, or visit us at www.miltex.com.



Illustrations and content provide general description only and may be subject to change. Miltex, Miltex Redefining Excellence, N-Tralig and GripRite are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra LifeSciences Corporation. Copyright © 2009 Integra LifeSciences Corporation