

MILTEX HI-LIGHT

THOMPSON JUNIOR



The Miltex Thompson Jr. is designed to handle calculus in the most difficult to reach areas. The finely tipped scaler adapts to give you access interproximally in the posterior part of the mouth as well as in the anterior portion. It is an excellent instrument to use on patients with healthy tight tissue but who have slight to moderate subgingival calculus. The overall length of the tip allows for easy access for this often hard to reach calculus. The Thompson Jr. was the first instrument of its kind and is now available in GripLite, GripLite S6 and stainless steel size 6 handles.

Part #	Description
80-5008	Scaler, Thompson Jr., GripLite S6
70-5008	Scaler, Thompson Jr., GripLite
70-721	Scaler, Thompson Jr. Stainless Steel #6 Handle
TPTHJR	Scaler, Thompson Jr. Tactile Tone Handle

FEATURES

- Finely tipped scaler
- Longer tip
- Available in multiple handle configurations

BENEFITS

- Adapts interproximally in the posterior part of the mouth
- Provides access for subgingival calculus
- Meet the needs of dental professionals

Miltex
Redefining Excellence

phone 717 840-9335
toll free 866 854-8300
www.miltex.com

1007/0/A/1

DENT-THMJNR-HI



+H834DENTTHMJNRHION