

Miltex[®]
Redefining Excellence



DENTAL SURGICAL INSTRUMENTS
...the most popular patterns!

...for GENERAL INSTRUMENTS

SCISSORS

Our razor sharp scissors are available in a variety of ergonomic designs in both ring handle and Castroviejo patterns. From cutting delicate tissue to removing sutures to performing precision procedures, Miltex is recognized as a leader in delivering consistent cutting performance and durability.

NEEDLE HOLDERS

A close examination of Miltex® Needle Holders, either standard or tungsten carbide, will reveal meticulously crafted surgical tools that will perform precisely as intended. Smoothly beveled box locks and jaw edges protect delicate suture material from snags or cutting. The jaw surface is milled to ensure a secure grasp and to prevent needle slippage or rotation.

HEMOSTATS

Our hemostats are designed for use in a variety of clinical dental procedures from the routine to the most delicate. Miltex is dedicated to providing customers with a wide variety of quality instruments with just the right feel and performance. Available in a variety of styles and manufactured from the highest quality stainless steel to exacting specifications, our hemostats demonstrate that Miltex is dedicated to providing you with just the right instrument.

TISSUE & DRESSING FORCEPS

Miltex® tissue and dressing forceps are designed to handle the most delicate procedures as well as routine surgical tasks. We offer a comprehensive line of versatile tissue and dressing forceps to handle a variety of needs.



SCISSORS

IRIS

- All purpose scissors
- Straight and curved patterns available
- SuperCut & tungsten carbide patterns available



5-306

KELLY

- Available with one smooth blade and one serrated blade
- All purpose scissors
- Straight and curved patterns available



5-252

STITCH

- Straight and curved patterns available
- Wide range of lengths available
- Notched blade on the Spencer, Stitch and Littauer scissors helps lift suture for easy cutting



9-102

REF	Description
5-300	Iris Scissors 4 Straight
5-302	Iris Scissors 4 Curved
5D-303	Iris Operating Scissors 4-1/2 Curved sideways
5-304	Iris Scissors 4-1/2 Straight
5-304TC*	Iris Scissors 4-1/2 Straight, Tungsten Carbide
5-SC-304**	Iris Scissors 4-1/2 Straight, SuperCut
5-306	Iris Scissors 4-1/2 Curved
5-306TC*	Iris Scissors 4-1/2 Curved, Tungsten Carbide
5-SC-306**	Iris Scissors 4-1/2 Curved, SuperCut
5-308	Iris Scissors 4-1/2 Angled on side
18-1396	Iris Scissors 3-1/2 Straight, with 20 mm blades, delicate
18-1398	Iris Scissors 3-1/2 Curved, with 20 mm blades, delicate
5-250	KELLY Scissors 6-1/4 Straight, sharp points
5D-251	KELLY Scissors 6-1/4 Straight, sharp points, one serrated blade
5-252	KELLY Scissors 6-1/4 Curved, sharp points
5D-253	KELLY Scissors 6-1/4 Sharp points, curved, one serrated blade
5-254	KELLY Scissors 6-1/4 Sharp points, angled on side
5D-255	KELLY Scissors 6-1/4 Sharp points, angular, one serrated blade
9-100	Spencer Stitch Scissors 3-1/2
9-108	Stitch Scissors 4-1/2 Angled, Delicate
9-102	Littauer Stitch Scissors 4-1/2
9-104	Littauer Stitch Scissors 5-1/2
9-110	O'BRIEN Stitch Scissors 3-3/4 Angled, sharp points
9-115	REEH Stitch Scissors 3-7/8 sharp tips, small hook, one serrated blade
9-116	MILTEX® Stitch Scissors 3-7/8 Curved, sharp pointed tips
9-117	MILTEX® Stitch Scissors 4-7/8 Curved, sharp pointed tips
5-264	DEAN Dissecting Scissors 6-3/4 Angled on flat, one serrated blade
5-320	GOLDMAN-FOX Scissors Curved fine tips, one serrated blade
5-320TC*	GOLDMAN-FOX Scissors Curved fine tips, one serrated blade, tungsten carbide
5D-322	GOLDMAN-FOX Scissors Curved, one serrated blade, angled shank
5D-324	LA GRANGE Scissors 4-1/4 Curved, Angled, one serrated blade

IRIS

KELLY

STITCH

OTHER

5-SC-306**



5D-251



9-108



5-320TC*



5D-324



*Premium all purpose scissors with tungsten carbide inserts

**Premium tissue only scissors with razor-sharp upper blade and serrated lower blade

SCISSORS

CROWN & COLLAR

REF	Description
9D-115	RA Matrix Scissors 4-1/4 one serrated blade
9D-118	Pedo Crown-Festooning Scissors 3-1/2 Curved
9-120	Wire Cutting Scissors 4-3/4 Straight, one serrated blade
9D-121	Wire Cutting Scissors 4-3/4 Straight
9D-132	Wire Cutting Scissors 4-3/4 Straight, one serrated blade
9D-133	Wire Cutting Scissors 4-3/4 Straight
9D-134	Wire Cutting Scissors 4-1/4 Curved, one serrated blade
9D-135	Wire Cutting Scissors 4-1/4 Curved



9D-135

NEEDLE HOLDERS

STANDARD (STD)


Miltex® needle holders are made from the finest quality German stainless steel, which ensures a high quality product.



8-40

TUNGSTEN CARBIDE (TC)

Tungsten carbide needle holders offer longer wear and superior grasping power to prevent needle slippage.



8-30TC

STD	TC	Description
8-40	8-40TC	Baumgartner Needle Holder 5-1/2
8-8	8-8TC	Halsey Needle Holder 5 (smooth)
	8-8A-TC	Halsey Needle Holder 5 (serrated)
8-4	8-4TC	Derf Needle Holder 4-3/4
8-50	8-50TC	Crile-Wood Needle Holder 6
8-3	8-3TC	PAR Needle Holder 4-1/2, delicate serrated jaws
8-5	8-5TC	PAR Needle Holder 5, delicate serrated jaws
8-28	8-28TC	MATHIEU Needle Holder 5-1/2, spring handles
8-30	8-30TC	MATHIEU Needle Holder 6-1/4, spring handles
8-42	8-42TC	Mayo-Hegar Needle Holder 5
8-44	8-44TC	Mayo-Hegar Needle Holder 6
8-46	8-46TC	Mayo-Hegar Needle Holder 7
8-15	8-15TC	Olsen-Hegar Needle Holder 5-1/2 (with suture scissors)
8-16	8-16TC	Olsen-Hegar Needle Holder 6-1/2 (with suture scissors)

*Miltex can replace tungsten carbide insert

8-40



8-8



8-4



8-50



8-5TC



8-42TC



8-16



OLSEN-HEGAR Needle Holders are available in a combination pattern with suture scissors.

HEMOSTATS

STANDARD (STD)



7-4

HEMOSTATS

STD	Description
7-4	Halsted Mosquito Forceps 5 Curved
7-2	Halsted Mosquito Forceps 5 Straight
7-20	Halsted Mosquito Forceps 5 Angled 45° on flat
7-26	Hartman Mosquito Forceps 3-1/2 Curved
7-24	Hartman Mosquito Forceps 3-1/2 Straight
7-38	Kelly Forceps 5-1/2 Curved
7-36	Kelly Forceps 5-1/2 Straight



TISSUE & DRESSING FORCEPS


STANDARD (STD)



6-118

TISSUE & DRESSING FORCEPS

STD	Description
6-118	Adson Dressing Forceps 4-3/4 (serrated)
6-120	Adson Tissue Forceps 4-3/4 (1x2)
6-121	Adson Tissue Forceps 4-3/4 (1x2 and cross serrated)
6-123	Adson Tissue & Suture Forceps 4-3/4 (1x2 with tying platform)
6-124	Brown-Adson Tissue Forceps 4-3/4 (7x7 side grasping)
6-112	Hudson Dressing Forceps 4-3/4 (serrated)
6-114	Hudson Tissue Forceps 4-3/4 (1x2)
6-4	Dressing Forceps 4-1/2 (serrated)
6-6	Dressing Forceps 5 (serrated)
6-8	Dressing Forceps 5-1/2 (serrated)
6-10	Dressing Forceps 6 (serrated)
6-40	Tissue Forceps 4-1/2 (1x2)
6-42	Tissue Forceps 5 (1x2)
6-44	Tissue Forceps 5-1/2 (1x2)
6-104	Semkin Dressing Forceps 5 (serrated)



Additional patterns are available...
Contact your local Miltex Representative



www.miltex.com
toll-free phone 866-854-8300
toll-free fax 866-854-8400
phone 717-840-9335
fax 717-840-9347

Please call Customer Service today for additional information, an in-office demonstration, to place an order, or to receive other Miltex catalogs.

Miltex and Miltex Redefining Excellence are registered trademarks of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. The Integra wave logo is a trademark of Integra LifeSciences Corporation. Copyright© 2010 Integra LifeSciences Corporation

