


MILTEX HI-LIGHT

MICRO BOLEY GAUGE



REF	Description
68-697	Micro Boley Gauge, Stainless Steel

The Miltex Micro Boley Gauge is used intraorally to measure tooth dimensions and interproximal space for cosmetic/restorative and orthodontic procedures. This sharp-looking instrument features slender beaks to reach small spaces in the mouth and can make accurate inside dimension (ID) or outside dimension (OD) measurements. The gauge can measure up to 95mm while the measurement markings are easily readable. A thumb screw locks the sliding beak to ensure accurate readings. Additionally, the Micro Boley Gauge also features a tip that measures tooth depth for crown and bridge restoratives.

Miltex also offers other gauges and bow dividers used for cosmetic and restorative procedures.

REF	Description
68-694	Boley Gauge, Stainless Steel
68-693	Crown Gauge, with 2 Sided Scale
68-695	Spring Bow Divider, 4.5" Radius
68-696	Spring Bow Divider, 6" Radius

Miltex
Redefining Excellence

phone 717 840-9335
toll free 866 854-8300
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.