

# MILTEX HI-LIGHT

## PETITE SYRINGE



REF	Description
77-70	Petite Aspirating Syringe, Cook-Waite (CW)

The Miltex® Petite Syringe was especially designed to meet the needs of dental professionals who prefer a more comfortable, form-fitting syringe. Ergonomical by design, the Miltex Petite Syringe features:

- **Reduced thumb ring handle diameter** to eliminate the excess gap, allowing just enough room for comfort and ease.
- **New oval shaped ring** handle that fits naturally to the shape of your thumb.
- **Rounded edges** on ring handle to create a smooth and soft feel.
- **Slightly flared base** on ring handle to ease pressure while dispensing.
- **Shortened harpoon rod** to reduce strain during aspiration.
- **Wide wings** allow for a better grip with greater control.

*Miltex*  
Redefining Excellence

phone 717 840-9335  
toll free 866 854-8300  
www.miltex.com

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email [customerservice@miltex.com](mailto:customerservice@miltex.com), or visit us at [www.miltex.com](http://www.miltex.com).

