

MILTEX® HI-LIGHT

SHORT HANDLED BIOPSY PUNCH



The Short Handled Biopsy Punch is designed to enable dentists to take tissue samples in all areas of the mouth. It features the same razor-sharp blade as the original Miltex® biopsy punches, but the shorter handle was modified to allow dental professionals access in even the hardest to reach, posterior areas of the mouth. The total length of the instrument is a mere 1-1/16 inch!

The product is provided in sterile, individually wrapped packets for easy use and is available in the four most common sizes: 3.5mm, 4mm, 5mm and 6mm.

REF	Description	Quantity
33-33-SH	Short Handled Biopsy Punch, 3.5mm	50/box
33-34-SH	Short Handled Biopsy Punch, 4mm	50/box
33-35-SH	Short Handled Biopsy Punch, 5mm	50/box
33-36-SH	Short Handled Biopsy Punch, 6mm	50/box



3.5mm



4mm



5mm



6mm

Illustrations and content signify general description only and may be subject to minor changes. For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.

Miltex is a registered trademark of Miltex, Inc., a subsidiary of Integra LifeSciences Corporation. Miltex Redefining Excellence is a trademarks of Miltex, Inc. The Integra Wave Logo is a trademark of Integra Life Sciences Corporation. Copyright © 2008 Integra LifeSciences Corporation



+H834DENTSHBPHI0A