


Miltex[®]
Redefining Excellence



MILTEX DENTAL MOUTH MIRRORS

Miltex

Miltex 80-5450
Stainless USA Lot no. K3 CE

DenLite® Illuminated Mirror

The DenLite Illuminated Mirror System is designed to give Dental Professionals a shadow-free view of the mouth during routine dental exams. It provides Dental Professionals with a bright white light, which eliminates the need for constant adjustment of cumbersome overhead lights during everyday procedures. The DenLite system is the only cordless, illuminated dental instrument that can be used as a mirror while also being strong enough to be also used as a retractor during an examination.

In May of 2006 the DenLite Illuminated Mirror System, designed by Welch Allyn, won the prestigious Medical Design Excellence Award. This award is presented by Cannon Communications LLC (Los Angeles) and is the only awards program that exclusively recognizes contributions and advancements in the design of medical products.



See The Light...

Bright White Light - a one-watt LED light provides shadow-free, glare-free illumination that's safe for use around composites

Convenience - a lightweight cordless ergonomic handpiece

Versatility - split-beam technology and rugged design make it ideal for use as a mirror and as a retractor

Safety - disposable mirror heads and autoclavable hand piece eliminate the risk of cross-contamination

