

# MILTEX HI-LIGHT

## HI-5 ENDODONTIC FILES



Calcified root canals can be the most challenging and difficult endodontic procedures. Traditional K-files can be too flexible to break through calcification, but Miltex Hi-5 files are designed with enhanced tensile strength, allowing easier location, penetration and initial instrumentation of even the most complex or calcified canals.

Hi-5 files have five helical ground flutes, creating a pentagonal cross-section for maximum strength and penetration without bending or folding. Hi-5 files have a shorter pitch with more spirals than a conventional instrument, producing a very stiff, resilient and effective cutting instrument.

### FEATURES

- Penetrates calcified canals
- Unique pentagon shape is more stiff compared to K-files
- Short pitch helical flutes enhance stiffness and cutting efficiency
- Packaged 6 files per pack

REF	Description (6/pk)	Size
012-26001	Hi-5 Penetrating Hand File	06 21 mm
012-26002	Hi-5 Penetrating Hand File	08 21 mm
012-26003	Hi-5 Penetrating Hand File	10 21 mm
012-26004	Hi-5 Penetrating Hand File	15 21 mm
012-26005	Hi-5 Penetrating Hand File	06 25 mm
012-26006	Hi-5 Penetrating Hand File	08 25 mm
012-26007	Hi-5 Penetrating Hand File	10 25 mm
012-26008	Hi-5 Penetrating Hand File	15 25 mm
012-26009	Hi-5 Penetrating Hand File	Assorted 21 mm*
012-26010	Hi-5 Penetrating Hand File	Assorted 25 mm*

\*Assorted pack shown in photo above.

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email [customerservice@miltex.com](mailto:customerservice@miltex.com), or visit us at [www.miltex.com](http://www.miltex.com).

*Miltex*  
Redefining Excellence

toll free 866 854-8300  
fax 866 854-8400  
[www.miltex.com](http://www.miltex.com)