

# MILTEX HI-LIGHT

## MARK VII APEX LOCATOR

RELIABLE, COMPACT, ACCURATE



The Mark VII Apex Locator provides precision measurement of the canal length under all conditions including dry, wet and bleeding. Precisely measuring the distance from the tip of the endodontic file to the foramen. The Mark VII allows dental professionals to know the depth of the canal by showing when they are 2.0 mm, 1.5 mm, 1.0 mm, 0.5 mm, 0.25 mm away from the apex.

The Mark VII Apex Locator is easy to use and designed to be placed near the oral cavity, in the operator's field of vision. The apex locator is often clipped to the patient's napkin to ensure that the unit does not accidentally fall. The Mark VII does not require calibration between measuring multiple canals, thereby immediately enabling accurate measurement of another canal.

Part Number	Description
017-23180	Mark VII Apex Locator

*Miltex*  
Redefining Excellence

phone 717 840-9335  
toll free 866 854-8300  
www.miltex.com

0907/10/B/2

DENT-MKVII-HI



\*+H834DENTMKVIIHI0 \*