

MILTEX HI-LIGHT

XCISION® FORCEPS



The Miltex® Xcision® brand of forceps are the original forceps featuring a unique, light weight, ergonomic handle. The circular openings help to reduce hand and wrist fatigue and make the instrument easier to clean. The Xcision® forceps have a wear resistant Teflon disc which helps reduce wear in the joint of the instruments. They also feature a slim, compact design for easy access to molars.

The Miltex® Xcision® extracting forceps have been designed to allow subgingival access to the root while removing a tooth. The beaks are anatomically shaped and tapered, which allows them to penetrate deeper into the tooth socket and maximize the tooth-to-forcep contact. Serrated beaks allow for better gripping of the tooth surface, which helps aid in the atraumatic extraction of the tooth and reduces the risk of fracturing the root structure or causing trauma to the surrounding hard and soft tissues.

Features	Benefits
Tapered beaks	Allow access subgingivally into the socket of the tooth Aid in atraumatic removal of tooth
Anatomically shaped beaks	Perfectly match teeth
Serrated beaks	Better gripping of tooth surface Aid in atraumatic removal of tooth
High quality German Stainless Steel	Product wears well and will last a long time, maximizing instrument investment
Smooth handle and circular openings	Easier to clean vs. serrated handles
Lightweight stainless steel instrument	Reduces hand and wrist fatigue
Teflon disc	Helps reduce wear in the joint of the instrument

REF	Description
DEFXC7	Extracting Forceps, Xcision® 7, Upper Premolar Universal
DEFXC17	Extracting Forceps, Xcision® 17, Upper Molars, Right
DEFXC18	Extracting Forceps, Xcision® 18, Upper Molars, Left
DEFXC321	Extracting Forceps, Xcision® 321, Lower Premolars
DEFXC34	Extracting Forceps, Xcision® 34, Upper Anteriors
DEFXC35	Extracting Forceps, Xcision® 35, Upper Premolars
DEFXC35M	Extracting Forceps, Xcision® 35M, Upper Premolars, Modified Beak
DEFXC36	Extracting Forceps, Xcision® 36, Lower Anteriors
DEFXC36M	Extracting Forceps, Xcision® 36M, Lower Anteriors, Modified Beak
DEFXC44	Extracting Forceps, Xcision® 44, Upper Roots Universal
DEFXC45	Extracting Forceps, Xcision® 45, Lower roots
DEFXC45S	Extracting Forceps, Xcision® 45S, Lower Roots with Serrations
DEFXC51SB	Extracting Forceps, Xcision® 51SB, Upper Roots Universal
DEFXC51S	Extracting Forceps, Xcision® 51S, Upper Roots, With Serrated Edges
DEFXC51W	Extracting Forceps, Xcision® 51W, Upper Roots, With Serrated Beaks
DEFXC67	Extracting Forceps, Xcision® 67, Upper 3rd Molars
DEFXC67L	Extracting Forceps, Xcision® 67L, Upper 3rd Molars Left
DEFXC67R	Extracting Forceps, Xcision® 67R, Upper 3rd Molars Right
DEFXC73	Extracting Forceps, Xcision® 73, Lower Molars
DEFXC79	Extracting Forceps, Xcision® 79, Lower 3rd Molars
DEFXC79M	Extracting Forceps, Xcision® 79M, Lower 3rd Molars Modified Beak
DEFXC23	Extracting Forceps, Xcision® 23, Cowhorn Lower Molars
DEFXC3S	Extracting Forceps, Xcision® Crown Spreading

See reverse side for tip details

For additional information, an in-office demonstration, or to place an order, please contact your local Miltex, Inc. authorized Distributor or Miltex, Inc. Customer Service at 1-866-854-8300, email customerservice@miltex.com, or visit us at www.miltex.com.





DEFXC7 Extracting Forceps, Xcision® 7, Upper Premolar Universal



DEFXC17 Extracting Forceps, Xcision® 17, Upper Molars, Right



DEFXC18 Extracting Forceps, Xcision® 18, Upper Molars, Left



DEFXC321 Extracting Forceps, Xcision® 321 Lower Premolars



DEFXC34 Extracting Forceps, Xcision® 34, Upper Anteriors



DEFXC35 Extracting Forceps, Xcision® 35, Upper Premolars



DEFXC35M Extracting Forceps, Xcision® 35M, Upper Premolars, Modified Beak



DEFXC23 Extracting Forceps, Xcision® 23, Cowhorn Lower Molars



DEFXC36 Extracting Forceps, Xcision® 36, Lower Anteriors



DEFXC36M Extracting Forceps, Xcision® 36M, Lower Anteriors, Modified Beak



DEFXC44 Extracting Forceps, Xcision® 44, Upper Roots Universal



DEFXC45 Extracting Forceps, Xcision® 45, Lower roots



DEFXC45S Extracting Forceps, Xcision® 45S, Lower Roots with Serrations



DEFXC51S Extracting Forceps, Xcision® 51S, Upper Roots, With Serrated Edges



DEFXC51W Extracting Forceps, Xcision® 51W, Upper Roots, With Serrated Beaks



DEFXC51SB Extracting Forceps, Xcision® 51SB, Upper Roots Universal



DEFXC67 Extracting Forceps, Xcision® 67, Upper 3rd Molars



DEFXC67L Extracting Forceps, Xcision® 67L, Upper 3rd Molars Left



DEFXC67R Extracting Forceps, Xcision® 67R, Upper 3rd Molars Right



DEFXC79 Extracting Forceps, Xcision® 79, Lower 3rd Molars



DEFXC79M Extracting Forceps, Xcision® 79M, Lower 3rd Molars Modified Beak



DEFXC73 Extracting Forceps, Xcision® 73, Lower Molars



DEFXCCS Extracting Forceps, Xcision® Crown Spreading