

## PerioStar® 3000 Archbow Sharpening Settings

ARCHBOW	INSTRUMENT	#
<b>Standard Setting</b>		
1	GRACEY CURETTE	1/2
1	GRACEY CURETTE	3/4
1	GRACEY CURETTE	5/6
2	GRACEY CURETTE	7/8
1	GRACEY CURETTE	9/10
1	GRACEY CURETTE	11/12
1	GRACEY CURETTE	13/14
1	GRACEY CURETTE	11/13
1	GRACEY CURETTE	12/14
1	GRACEY CURETTE	15/16
1	GRACEY CURETTE	17/18
1	LANGER CURETTE	1/2
1	LANGER CURETTE	3/4
1	LANGER CURETTE	5/6
1	LANGER CURETTE	17/18
1	COLUMBIA CURETTE	2L/2R
1	COLUMBIA CURETTE	4L/4R
1	COLUMBIA CURETTE	13/14
1	RULE	
3	SCALER	3/4
3	SCALER	SM00
3	SCALER	204
3	SCALER	204S
3	SCALER	204SD
3	SCALER	H204S
4	SCALER	H6/H7
4	SCALER	U15/33
4	SCALER	SH6/7
4	SCALER	H5/33